

MECP2 (phospho Ser80) antibody

Cat. No. GTX82658

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Human, Mouse, Rat

[References \(2 \)](#)[Package](#)

100 µl

Applications

Application Note

Optimal working dilutions should be determined experimentally by the end user.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM HEPES, 150mM NaCl, 0.01% BSA, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser80 conjugated to KLH
Purification	Serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx82658 WB Image

Western blot of neonatal rat brain showing specific immunolabeling of the ~75k MeCP2 protein phosphorylated at Ser 80 using MeCP2 (phospho Ser80) antibody (GTx82658).



For full product information, images and publications, please visit our [website](#).