

IFNAR1 (phospho Ser535/Ser539) antibody

Cat. No. GTX82662

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Rat

Package
100 µl

Applications

Application Note

Optimal working dilutions should be determined experimentally by the end user.

Calculated MW 64 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM HEPES, 150mM NaCl, 0.01% BSA, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser535/539 conjugated to KLH
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX82662 WB Image

Western blot of immunoprecipitates from HEK 293 cells transfected with

1. Mock
2. IFNAR1 WT and
3. IFNAR1 Ser 535 Arg and Ser 539 Arg mutants.

The immunolabeling is absent in the IFNAR1 Ser 535 and Ser 539 mutants. The labeling is blocked by the phosphopeptide (Phos) used as the antigen but not by the dephosphopeptide (Dephos).



For full product information, images and publications, please visit our [website](https://www.genetex.com).