

GABA A Receptor gamma 2 antibody

Cat. No. GTX82677

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Rat

Package

50 µl

Applications

Application Note

WB:1:1000

IHC:1:400.

Optimal working dilutions should be determined experimentally by the end user.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide corresponding to amino acid residues specific to the gamma 2 subunit conjugated to KLH
Purification	Neat Serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX82677 WB Image

Western blot of 10 ug of rat hippocampal lysate showing specific immunolabeling of the ~46k gamma 2-subunit of the GABA_A-R using GABA_A Receptor gamma 2 antibody (GTX82677).



For full product information, images and publications, please visit our [website](#).