

Rab25 antibody [3F12]

Cat. No. GTX82785

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human, Mouse

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of Rab25 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX82785 WB Image

WB analysis of MCF-7 (1), T47D (2) and GC7901 (3) cell lysate using GTX82785 Rab25 antibody [3F12].




GTX82785 ICC/IF Image

ICC/IF analysis of A549 cells using GTX82785 Rab25 antibody [3F12].

Green : Rab25


Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments



GTX82785 IHC-P Image

IHC-P analysis of esophagus tissue (left) and human lung cancer (right) using GTX82785 Rab25 antibody [3F12].



GTX82785 FCM Image

FACS analysis of NIH3T3 cells using GTX82785 Rab25 antibody [3F12].

Green : Rab25

Purple : negative control



For full product information, images and publications, please visit our [website](https://www.genetex.com).