

PTH1R antibody [4D2]

Cat. No. GTX82794

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 66 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human PTH1R expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX82794 ICC/IF Image

ICC/IF analysis of SK-BR-3 cells using GTX82794 PTH1R antibody [4D2].


Green : PTH1R

Blue: DRAQ5 fluorescent DNA dye.



GTX82794 WB Image

WB analysis of PTH1R (AA: 27-188)-hlgGfc transfected HEK293 cell lysate using GTX82794 PTH1R antibody [4D2].



GTX82794 IHC-P Image

IHC-P analysis of human kidney tissue using GTX82794 PTH1R antibody [4D2].



For full product information, images and publications, please visit our [website](https://www.genetex.com).