

IP3 Receptor I antibody

Cat. No. GTX82887

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-Fr, FCM, IP, IHC
Reactivity	Human, Mouse, Rat, Dog, Hamster, Pig

[References \(2 \)](#)
[Package](#)

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1,000
ICC/IF	2µg/ml
IHC-Fr	1:250
FCM	3-5 µg/10 ⁶ cells
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 313 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.3 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic Peptide: N(1829) K K K D D E V D R D A P S R K K A K E(1848)
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX82887 IHC-Fr Image**

IHC-Fr analysis of rat cerebellum tissue using GTX82887 IP3 Receptor I antibody.

**GTX82887 ICC/IF Image**

ICC/IF analysis of SH-SY5Y cells using GTX82887 IP3 Receptor I antibody. Panel e is a no primary antibody control.

Green : Primary antibody

Blue : Nuclei

Red : Actin

Fixation : 4% paraformaldehyde

Permeabilization : 0.1% Triton X-100 for 10 minute

Dilution : 2 µg/ml in 0.1% BSA and incubated for 3 hours at room temperature



For full product information, images and publications, please visit our [website](#).