

ECM29 antibody

Cat. No. GTX82906

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC
Reactivity	Human, Mouse

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:250
ICC/IF	1:250
IHC	Assay dependent

Not tested in other applications.

Calculated MW 204 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide corresponding to the C-terminal residues R(1821) P E L Q E K A A L L K K T L E N L E(1839) of human Ecm29.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX82906 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX82906 ECM29 antibody. Panel e is a no primary antibody control.

Green : Primary antibody


Blue : Nuclei

Red : Actin

Fixation : 4% paraformaldehyde

Permeabilization : 0.1% Triton X-100 for 10 minute


Dilution : 1:250 dilution in 0.1% BSA and incubated for 3 hours at room temperature



GTX82906 WB Image

WB analysis of whole cell extracts (30 µg lysate) of HeLa (Lane 1), K-562 (Lane 2) and NTERA-2 cl.D1 (Lane 3) using GTX82906 ECM29 antibody.

Dilution : 1:250



GTX82906 WB Image

WB analysis of HeLa cells using GTX82906 ECM29 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).