

Cytokeratin 19 antibody [CK 202]

Cat. No. GTX82955

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, IHC
Reactivity	Human

Package
250 µl

Applications

Application Note

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears) For IHC dilute conc. antibodies 1:50-1:100, use streptavidin~biotin system or polymer system, incubate 30 minutes at room temperature.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human mammary epithelium, solubilized by detergent.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).