

Vimentin antibody

Cat. No. GTX82958

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	IHC
Reactivity	Human

Package
500 µl

APPLICATION

Application Note

Immunohistochemistry (IHC) 1:50-1:100 using polymer system or streptavidin-biotin system. Formalin-fixed paraffin-embedded (FFPE) tissue sections and cell smears do not require antigen retrieval.

Calculated MW 54 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	10mM PO ₄ , 200mM NaCl 1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	3 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human Vimentin
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).