

Insulin antibody (Ready-to-use)

Cat. No. GTX82959

Host	Guinea pig
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
7 ml

Applications

Application Note

For IHC-P: Supplied in ready to use form. Optimal dilutions/concentrations should be determined by the end user.

Properties

Form	Liquid
Buffer	Total immunoglobulin fraction in universal antibody dilution buffer.
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Purified human insulin conjugated to carrier protein
Purification	Purified
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).