

## Chicken Lysozyme protein (Biotin)

**Cat. No. GTX82960-pro**

<b>Applications</b>	Functional Assay
<b>Species</b>	Chicken

References ( 6 )  
Package  
1 ml

## Applications

**Application Note**

Specific activity: > 20 units/mg.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	10mM PO <sub>4</sub> , 150mM NaCl
<b>Preservative</b>	0.05% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Expression System</b>	Egg whites
<b>Conjugation</b>	Biotin
<b>Note</b>	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).