

Chicken Lysozyme protein (Biotin)

Cat. No. GTX82960-pro

Applications Functional Assay

Species Chicken

References (5)

Package

1 ml

Applications

Application Note

Specific activity: > 20 units/mg.

Properties

Form Liquid

Buffer 10mM PO₄, 150mM NaCl

Preservative 0.05% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 5 mg/ml (Please refer to the vial label for the specific concentration.)

Expression System Egg whites

Conjugation Biotin

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).