

Chicken Lysozyme protein (Biotin)

Cat. No. GTX82960-pro

Application	Functional Assay
Species	Chicken

Reference (3)
 Package
 1 ml

APPLICATION

Application Note

Specific activity: > 20 units/mg.

PROPERTIES

Form	Liquid
Buffer	10mM PO ₄ (pH6.2), 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5 mg/ml (Please refer to the vial label for the specific concentration.)
Expression System	Egg whites
Conjugation	Biotin
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).