

## Sialosyl-Tn antigen antibody [STn 219]

Cat. No. GTX82967

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, IHC-P, IHC
Reactivity	Human

References ( 1 )

Package

500 µl

## Applications

## Application Note

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears. For IHC dilute conc. antibodies 1:50-1:100, use streptavidin~biotin system or polymer system, incubate 30 minutes at room temperature.

## Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Ovine submaxillary mucin (OSM)
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).