

Tn antigen antibody [Tn 218]

Cat. No. GTX82968

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	ICC/IF, IHC-P, IHC-Fr
Reactivity	Human, Mouse

References (3)

Package

500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	1:50-1:100
IHC-Fr	1:50-1:100

Not tested in other applications.

Product Note

Epitope: alpha GalNacOSer/Thr

This antibody stains Tn antigen in epithelial carcinoma, pancreatic, colon, stomach, and breast.

Properties

Form	Liquid
Buffer	tissue culture medium
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Asialo ovine submaxillary mucin (AOSM).
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).