

## CD45 antibody [IH1]

**Cat. No. GTX82980**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Applications</b>	ICC/IF, IHC-Fr, IP
<b>Reactivity</b>	Human

**Package**  
500 µg

## Applications

**Application Note**

For IF: Use at a concentration of 5-20µg/ml. For IP: Use at a concentration of 5-20µg/ml. Optimal dilutions/concentrations should be determined by the researcher.

**Product Note** Detects all isoforms (180-220 kDa)

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Leukaemic lymphocytes from guinea pig
<b>Purification</b>	Protein A affinity purified
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).