

CD45 antibody [IH1]

Cat. No. GTX82980

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Application	ICC/IF, IHC-Fr, IP
Reactivity	Human

Package
500 µg

APPLICATION

Application Note

For IF: Use at a concentration of 5-20µg/ml. For IP: Use at a concentration of 5-20µg/ml. Optimal dilutions/concentrations should be determined by the researcher.

Product Note Detects all isoforms (180-220 kDa)

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Leukaemic lymphocytes from guinea pig
Purification	Protein A affinity purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).