

## IRE1 alpha antibody [9F2]

Cat. No. GTX82996

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 110 kDa. ([Note](#))

## Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human IRE1 alpha (aa282-433) expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




**GTX82996 IHC-P Image**

IHC-P analysis of human Placenta tissue using GTX82996 IRE1 alpha antibody [9F2].



**GTX82996 IHC-P Image**

IHC-P analysis of human brain tissue (A) and stomach tissue (B) using GTX82996 IRE1 alpha antibody [9F2].



**GTX82996 WB Image**

WB analysis of Raji (1), A431 (2), Jurkat (3), HeLa(4) and HEK293 (5) cell lysate using GTX82996 IRE1 alpha antibody [9F2].



For full product information, images and publications, please visit our [website](https://www.genetex.com).