

CD18 antibody [10E12]

Cat. No. GTX83119

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ICC/IF, ELISA
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of CD18 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx83119 ICC/IF Image

ICC/IF analysis of HL60 cells using GTx83119 CD18 antibody [10E12].

Green : CD18

Blue: DRAQ5 fluorescent DNA dye



GTx83119 ICC/IF Image

ICC/IF analysis of BCBL-1 cells (left) and L1210 cells (right) using GTx83119 CD18 antibody [10E12].

Red : CD18

Blue: DRAQ5 fluorescent DNA dye



For full product information, images and publications, please visit our [website](#).