

FGFR4 antibody [7H1]

Cat. No. GTX83129

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 88 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant extracellular fragment of human FGFR4 fused with hlgGfc tag expressed in HEK293 cell line.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83129 ICC/IF Image

ICC/IF analysis of methanol-fixed HEK293 cells transfected with FGFR4-hlgFc using GTX83129 FGFR4 antibody [7H1].

Green : FGFR4

Blue: DRAQ5 fluorescent DNA dye


GTX83129 WB Image

WB analysis of extracellular domain of human FGFR4 (aa22-369) using GTX83129 FGFR4 antibody [7H1].



For full product information, images and publications, please visit our [website](https://www.genetex.com).