

Emerin antibody

Cat. No. GTX83131

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Mouse, Rat

Package
50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:1000

For ICC/IF: Use at 1:100

For EM: Use at 1:50

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	29 kDa. (Note)
Product Note	Rodent

Properties

Form	Liquid
Buffer	PBS, 10% Gelatin
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Rat emerin residues 114-183 (GST tag at Nterminus)
Purification	Affinity purified on GST-sepharose column with attached antigen
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).