

Emerin antibody

Cat. No. GTX83131

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Mouse, Rat

Package
50 µl

Applications

Application Note**Recommended Starting Dilutions:**

For WB: Use at 1:1000

For ICC/IF: Use at 1:100

For EM: Use at 1:50

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 29 kDa. ([Note](#))**Product Note** Rodent

Properties

Form Liquid**Buffer** PBS, 10% Gelatin**Preservative** 0.1% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Immunogen** Rat emerin residues 114-183 (GST tag at Nterminus)**Purification** Affinity purified on GST-sepharose column with attached antigen**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).