

## FRK antibody [1A8H2C12]

Cat. No. GTX83187

|                     |                  |
|---------------------|------------------|
| <b>Host</b>         | Mouse            |
| <b>Clonality</b>    | Monoclonal       |
| <b>Isotype</b>      | IgG1             |
| <b>Applications</b> | WB, IHC-P, ELISA |
| <b>Reactivity</b>   | Human            |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1/500 - 1/2000       |
| IHC-P              | 1/200 - 1/1000       |
| ELISA              | 1/10000              |

Not tested in other applications.

**Calculated MW** 58 kDa. ( [Note](#) )

## Properties

|                     |                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>         | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>       | Ascites                                                                                                                                                                                                                    |
| <b>Preservative</b> | 0.03% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>      | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Immunogen</b>    | Purified recombinant fragment of FRK (aa2-300) expressed in E. Coli.                                                                                                                                                       |
| <b>Purification</b> | Unpurified                                                                                                                                                                                                                 |
| <b>Conjugation</b>  | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX83187 WB Image

WB analysis of HEK293T cells transfected with the pCMV6-ENTRY control (1) and pCMV6-ENTRY FRK cDNA (2) using GTX83187 FRK antibody [1A8H2C12].



### GTX83187 IHC-P Image

IHC-P analysis of human breast tissue using GTX83187 FRK antibody [1A8H2C12].



For full product information, images and publications, please visit our [website](https://www.genetex.com).