

Cytokeratin 15 antibody [6E7]

Cat. No. GTX83193

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of KRT15 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83193 WB Image


WB analysis of A431 cell lysate using GTX83193 Cytokeratin 15 antibody [6E7].


GTX83193 FCM Image

FACS analysis of PACN-1 cells using GTX83193 Cytokeratin 15 antibody [6E7].

Green : Cytokeratin 15

Purple : negative control



GTX83193 ICC/IF Image

ICC/IF analysis of HepG2(left) and PACN-1 (right) cells using GTX83193 Cytokeratin 15 antibody [6E7].

Green : Cytokeratin 15

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments


GTX83193 IHC-P Image

IHC-P analysis of human Tonsil tissue using GTX83193 Cytokeratin 15 antibody [6E7].



For full product information, images and publications, please visit our [website](https://www.genetex.com).