

GABPA antibody [8C1B10]

Cat. No. GTX83195

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Human, Mouse

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IP	Assay dependent
ELISA	1/10000

Not tested in other applications.

Calculated MW 51 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human GABPA (aa120-190) expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX83195 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX83195 GABPA antibody [8C1B10].


Green : GABPA

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments


GTX83195 WB Image

WB analysis of HEK293T cells transfected with the pCMV6-ENTRY control (1) and pCMV6-ENTRY GABPA cDNA (2) using GTX83195 GABPA antibody [8C1B10].


GTX83195 WB Image

WB analysis of HeLa (1), A549 (2), MCF-7 (3), NIH3T3 (4) and SMMC-7721 (5) cell lysate using GTX83195 GABPA antibody [8C1B10].



For full product information, images and publications, please visit our [website](https://www.genetex.com).