

KSHV K8 alpha antibody [8C12G10G1]

Cat. No. GTX83237

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, ELISA
Reactivity	Kaposi's sarcoma-associated herpesvirus

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of KSHV K8α expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX83237 WB Image

WB analysis of BCBL-1 (1) and TPA induced BCBL-1 (2) cell lysate using GTX83237 KSHV K8 alpha antibody [8C12G10G1].



GTX83237 ICC/IF Image

ICC/IF analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) with AEC staining using GTX83237 KSHV K8 alpha antibody [8C12G10G1].



For full product information, images and publications, please visit our [website](https://www.genetex.com).