

LYRIC antibody [2F11C3]

Cat. No. GTX83259

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA, IHC
Reactivity	Human

References (2)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000
IHC	Assay dependent

Not tested in other applications.

Calculated MW 64 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human LYRIC(AA : 378-579) expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83259 FCM Image


FACS analysis of HeLa cells using GTX83259 LYRIC antibody [2F11C3].

Green : LYRIC

Purple : negative control


GTX83259 IHC-P Image

IHC-P analysis of human Liver tissue using GTX83259 LYRIC antibody [2F11C3].


GTX83259 WB Image

WB analysis of K562 (1), SKBR-3 (2), T47D (3), HeLa (4) and MCF-7 (5) cell lysate using GTX83259 LYRIC antibody [2F11C3].



For full product information, images and publications, please visit our [website](https://www.genetex.com).