

Androgen Receptor antibody [2H8]

Cat. No. GTX83309

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 99 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human AR(AA : 221-321) expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx83309 WB Image


WB analysis of Androgen receptor (aa221-321)-hlgGfc transfected HEK293 cell lysate using GTx83309 Androgen Receptor antibody [2H8].


GTx83309 IHC-P Image

IHC-P analysis of human Prostate tissue using GTx83309 Androgen Receptor antibody [2H8].


GTx83309 WB Image

WB analysis of K562 (1), Jurkat (2) and LNCaP (3) cell lysate using GTx83309 Androgen Receptor antibody [2H8].


GTx83309 IHC-P Image

IHC-P analysis of human human testis (left) and prostate cancer (right) tissue using GTx83309 Androgen Receptor antibody [2H8].



For full product information, images and publications, please visit our [website](https://www.genetex.com).