

nm23-H1 antibody [4B2]

Cat. No. GTX83325

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 17 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human NME1 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX83325 IHC-P Image

IHC-P analysis of human Placenta tissue using GTX83325 nm23-H1 antibody [4B2].




GTX83325 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX83325 nm23-H1 antibody [4B2].

Green : nm23-H1

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments




GTX83325 FCM Image

FACS analysis of Jurkat cells using GTX83325 nm23-H1 antibody [4B2].

Green : nm23-H1

Purple : negative control



GTX83325 WB Image

WB analysis of nm23-H1(AA: 27-177)-hlgGfC transfected HEK293 cell lysate using GTX83325 nm23-H1 antibody [4B2].



For full product information, images and publications, please visit our [website](https://www.genetex.com).