

Cathelicidin antibody

Cat. No. GTX83344

Host	Mouse
Clonality	Polyclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:1000

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	19 kDa. (Note)
Product Note	Human hCAP-18 115-NQARGSF-121 in pro domain

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Partially purified human neutrophil membrane proteins (heparin Ultrogel/lectin affinity purification fractions in octyl glucoside)
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).