

CD14 antibody [5A3]

Cat. No. GTX83395

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000


Not tested in other applications.

Properties


Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human CD14 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83395 IHC-P Image

IHC-P analysis of human lymphnode using GTX83395 CD14 antibody [5A3].


GTX83395 FCM Image

FACS analysis of human PBMC using GTX83395 CD14 antibody [5A3].



For full product information, images and publications, please visit our [website](#).