

TTLL12 antibody [12B7]

Cat. No. GTX83471

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:100
FCM	1:100

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.73 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human TTLL12 (NP_055955) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2


DATA IMAGES



GTX83471 WB Image

WB analysis of wild type and TTLL12 knock out 293T cell lysate(10µg per lane) using GTX83471 TTLL12 antibody [12B7].

Dilution : 1:500




GTX83471 FCM Image

FACS analysis of HeLa cells using GTX83471 TTLL12 antibody [12B7].

Red : Primary antibody


Blue : Negative control antibody



GTX83471 FCM Image

FACS analysis of HEK293T cells transfected with either TTLL12 plasmid (Red) or empty vector control plasmid (Blue) using GTX83471 TTLL12 antibody [12B7].

Dilution : 1:1000



GTX83471 WB Image

WB analysis of HEK293T cells transfected with TTLL12 plasmid (Right) or empty vector (Left) for 48 hrs using GTX83471 TTLL12 antibody [12B7].

Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).