

Melanoma gp100 antibody [7B2]

Cat. No. GTX83639

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100
FCM	1:100

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human SILV (NP_008859) produced in HEK293T cell.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83639 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with Melanoma gp100 plasmid using GTX83639 Melanoma gp100 antibody [7B2].


GTX83639 ICC/IF Image

ICC/IF analysis of HT29 cells using GTX83639 Melanoma gp100 antibody [7B2].


GTX83639 FCM Image

FACS analysis of HEK293T cells transfected with either Melanoma gp100 plasmid(Red) or empty vector control plasmid(Blue) using GTX83639 Melanoma gp100 antibody [7B2].

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).