

RND1 antibody [4F12]

Cat. No. GTX83709

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:100
FCM	1:100

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human RND1 (NP_055285) produced in HEK293T cell.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx83709 WB Image


WB analysis of HEK293T cells transfected with RND1 plasmid (Right) or empty vector (Left) for 48 hrs using GTx83709 RND1 antibody [4F12].

Loading : 5 ug per lane



GTx83709 FCM Image


FACS analysis of HEK293T cells transfected with either RND1 plasmid(Red) or empty vector control plasmid(Blue) using GTx83709 RND1 antibody [4F12].



GTx83709 ICC/IF Image

ICC/IF analysis of HeLa cell lysate using GTx83709 RND1 antibody [4F12].

Dilution : 1:200



GTx83709 WB Image

WB analysis of HeLa cell lysate using GTx83709 RND1 antibody [4F12].

Loading : 10 ug per lane

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).