

PSMD10 antibody [3F6]

Cat. No. GTX83773

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
FCM	1:100

Not tested in other applications.

Calculated MW 24 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human PSMD10 (NP_002805) produced in HEK293T cell.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX83773 WB Image


WB analysis of various cell lines using GTX83773 PSMD10 antibody [3F6].

Loading : 35 ug per lane


GTX83773 WB Image

WB analysis of HEK293T cells transfected with PSMD10 plasmid (Right) or empty vector (Left) for 48 hrs using GTX83773 PSMD10 antibody [3F6].


Loading : 5 ug per lane


GTX83773 FCM Image

FACS analysis of HeLa cells using GTX83773 PSMD10 antibody [3F6].

Red : Primary antibody

Blue : Negative control antibody


GTX83773 FCM Image

FACS analysis of Jurkat cells using GTX83773 PSMD10 antibody [3F6].

Red : Primary antibody

Blue : Negative control antibody



For full product information, images and publications, please visit our [website](https://www.genetex.com).