

# TrkC antibody [19B8]

**Cat. No. GTX83976**

|                     |                 |
|---------------------|-----------------|
| <b>Host</b>         | Mouse           |
| <b>Clonality</b>    | Monoclonal      |
| <b>Isotype</b>      | IgG1            |
| <b>Applications</b> | WB              |
| <b>Reactivity</b>   | Human, Rat, Dog |

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:2000               |

Not tested in other applications.

**Calculated MW** 94 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 1% BSA, 50% Glycerol                                                                                                                                                                                                  |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | Full-length protein expressed in 293T cell transfected with human NTRK3 expression vector                                                                                                                                  |
| <b>Purification</b>  | Purified by affinity chromatography                                                                                                                                                                                        |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


## DATA IMAGES



### GTx83976 WB Image

WB analysis of HEK293T cells transfected with TrkC plasmid (Right) or empty vector (Left) for 48 hrs using GTx83976 TrkC antibody [19B8].

Loading : 5 ug per lane



### GTx83976 WB Image

WB analysis of various cell lines using GTx83976 TrkC antibody [19B8].

Loading : 35 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).