

JNK1 antibody [1C2]

Cat. No. GTX84143

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-384 of human JNK1 (NP_002741) produced in E.coli.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx84143 WB Image

WB analysis of various cell lines using GTx84143 JNK1 antibody [1C2].

Loading : 35 ug per lane



GTx84143 WB Image


WB analysis of NIH/3T3 cell lysate using GTx84143 JNK1 antibody [1C2].

Dilution : 1:1000



GTx84143 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with JNK1 plasmid using GTx84143 JNK1 antibody [1C2].



GTx84143 WB Image

WB analysis of HEK293T cells transfected with JNK1 plasmid (Right) or empty vector (Left) for 48 hrs using GTx84143 JNK1 antibody [1C2].

Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).