

MEK1 antibody [1C1]

Cat. No. GTX84168

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:100

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full-length protein expressed in 293T cell transfected with human MAP2K1 expression vector
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX84168 ICC/IF Image


ICC/IF analysis of COS7 cells transiently transfected with MEK1 plasmid using GTX84168 MEK1 antibody [1C1].



GTX84168 WB Image

WB analysis of HEK293T cells transfected with MEK1 plasmid (Right) or empty vector (Left) for 48 hrs using GTX84168 MEK1 antibody [1C1].

Loading : 5 ug per lane



GTX84168 ICC/IF Image

ICC/IF analysis of HT29 cells using GTX84168 MEK1 antibody [1C1].



For full product information, images and publications, please visit our [website](#).