

LIN7B antibody [1C9]

Cat. No. GTX84216

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:100
IHC-P	1:150
FCM	1:100

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.58 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human LIN7B (NP_071448) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx84216 ICC/IF Image


ICC/IF analysis of COS7 cells transiently transfected with LIN7B plasmid using GTx84216 LIN7B antibody [1C9].


GTx84216 FCM Image

FACS analysis of HEK293T cells transfected with either LIN7B plasmid(Red) or empty vector control plasmid(Blue) using GTx84216 LIN7B antibody [1C9].


GTx84216 IHC-P Image

IHC-P analysis of human bladder carcinoma tissue using GTx84216 LIN7B antibody [1C9].
Antigen retrieval : Heat-induced epitope retrieval by 10mM citrate buffer, pH6.0, 100°C for 10min.
Dilution : 1:50


GTx84216 WB Image

WB analysis of HEK293T cells transfected with LIN7B plasmid (Right) or empty vector (Left) for 48 hrs using GTx84216 LIN7B antibody [1C9].
Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).