

GSTT2 antibody [2E3]

Cat. No. GTX84400

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	1:50-1:150
FCM	1:100

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.83 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human GSTT2 (NP_000845) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX84400 IHC-P Image


IHC-P analysis of human pancreas tissue using GTX84400 GSTT2 antibody [2E3].
Dilution : 1:150


GTX84400 IHC-P Image


IHC-P analysis of human breast adenocarcinoma tissue using GTX84400 GSTT2 antibody [2E3].
Dilution : 1:150


GTX84400 IHC-P Image

IHC-P analysis of human endometrium adenocarcinoma tissue using GTX84400 GSTT2 antibody [2E3].
Dilution : 1:150


GTX84400 IHC-P Image

IHC-P analysis of human liver tissue using GTX84400 GSTT2 antibody [2E3].
Dilution : 1:150


GTX84400 FCM Image

FACS analysis of HeLa cells using GTX84400 GSTT2 antibody [2E3].
Red : Primary antibody
Blue : Negative control antibody



For full product information, images and publications, please visit our [website](#).