

GSC antibody [4F11]

Cat. No. GTX84409

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.19 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment expressed in E.coli corresponding to amino acids 107-257 of human goosecoid
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX84409 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX84409 GSC antibody [4F11].



GTX84409 ICC/IF Image

ICC/IF analysis of HeLa cells transiently transfected with GSC plasmid using GTX84409 GSC antibody [4F11].



For full product information, images and publications, please visit our [website](#).