

FOXA1 antibody [3A8]

Cat. No. GTX84487

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF
Reactivity	Human, Dog

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-473 of human FOXA1
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX84487 WB Image

WB analysis of various cell lines using GTX84487 FOXA1 antibody [3A8].

Loading : 35 ug per lane


GTX84487 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with FOXA1 plasmid using GTX84487 FOXA1 antibody [3A8].


GTX84487 WB Image

WB analysis of HEK293T cells transfected with FOXA1 plasmid (Right) or empty vector (Left) for 48 hrs using GTX84487 FOXA1 antibody [3A8].

Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).