

## FKBP12 antibody [3B3]

Cat. No. GTX84492

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	WB, IHC-P, FACS
Reactivity	Human, Dog

Package  
100 µl

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	1:150
FACS	1:100

Not tested in other applications.

Calculated MW 12 kDa. ( [Note](#) )

## PROPERTIES

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.8 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human FKBP1A (NP_000792) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).




**DATA IMAGES**

**GTX84492 WB Image**


WB analysis of various cell lines using GTX84492 FKBP12 antibody [3B3].

Loading : 10 ug per lane

Dilution : 1:200



**GTX84492 FACS Image**

FACS analysis of HEK293T cells transfected with either FKBP12 plasmid(Red) or empty vector control plasmid(Blue) using GTX84492 FKBP12 antibody [3B3].


**GTX84492 IHC-P Image**

IHC-P analysis of human colon adenocarcinoma tissue using GTX84492 FKBP12 antibody [3B3].

Antigen retrieval : Heat-induced epitope retrieval by 10mM citrate buffer, pH6.0, 100°C for 10min.


**GTX84492 WB Image**

WB analysis of HEK293T cells transfected with FKBP12 plasmid (Right) or empty vector (Left) for 48 hrs using GTX84492 FKBP12 antibody [3B3].

Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).