

Fumarate hydratase antibody [2A2]

Cat. No. GTX84499

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human, Dog

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:50-1:100

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human FH (NP_000134) produced in HEK293T cell.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx84499 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with Fumarate hydratase plasmid using GTx84499 Fumarate hydratase antibody [2A2].


GTx84499 WB Image


WB analysis of HEK293T cells transfected with Fumarate hydratase plasmid (Right) or empty vector (Left) for 48 hrs using GTx84499 Fumarate hydratase antibody [2A2].

Loading : 5 ug per lane


GTx84499 WB Image

WB analysis of various cell lines using GTx84499 Fumarate hydratase antibody [2A2].

Loading : 35 ug per lane


GTx84499 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx84499 Fumarate hydratase antibody [2A2].



For full product information, images and publications, please visit our [website](https://www.genetex.com).