

EPM2AIP1 antibody [1F11]

Cat. No. GTX84561

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:100
IHC-P	1:50-1:150
FCM	1:100

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.61 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human EPM2AIP1 (NP_055620) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX84561 IHC-P Image

IHC-P analysis of human colon tissue using GTX84561 EPM2AIP1 antibody [1F11].
Dilution : 1:150




GTX84561 FCM Image

FACS analysis of HEK293T cells transfected with either EPM2AIP1 plasmid(Red) or empty vector control plasmid(Blue) using GTX84561 EPM2AIP1 antibody [1F11].



GTX84561 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with EPM2AIP1 plasmid using GTX84561 EPM2AIP1 antibody [1F11].



GTX84561 WB Image

WB analysis of HEK293T cells transfected with EPM2AIP1 plasmid (Right) or empty vector (Left) for 48 hrs using GTX84561 EPM2AIP1 antibody [1F11].
Loading : 5 ug per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).