

ICAD antibody [5C4]

Cat. No. GTX84621

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2500-1:5000
ICC/IF	1:50-1:100

Not tested in other applications.

Calculated MW 37 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human DFFA (NP_004392) produced in E.coli.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX84621 WB Image

WB analysis of HEK293T cells transfected with ICAD plasmid (Right) or empty vector (Left) for 48 hrs using GTX84621 ICAD antibody [5C4].

Loading : 5 ug per lane



GTX84621 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with ICAD plasmid using GTX84621 ICAD antibody [5C4].


GTX84621 WB Image

WB analysis of various cell lines using GTX84621 ICAD antibody [5C4].

Loading : 35 ug per lane


GTX84621 ICC/IF Image

ICC/IF analysis of A549 cells using GTX84621 ICAD antibody [5C4].



For full product information, images and publications, please visit our [website](https://www.genetex.com).