

Rho GDI alpha antibody [1A7]

Cat. No. GTX84860

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:100
IHC-P	1:50

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.72 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human ARHGDI (NP_004300) produced in HEK293T cell.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx84860 IHC-P Image

IHC-P analysis of kidney tissue using GTx84860 Rho GDI alpha antibody [1A7].

Antigen retrieval : Heat-induced epitope retrieval by 10mM citrate buffer, pH6.0, 100°C for 10min.

Dilution : 1:50




GTx84860 IHC-P Image

IHC-P analysis of colon adenocarcinoma tissue using GTx84860 Rho GDI alpha antibody [1A7].


Antigen retrieval : Heat-induced epitope retrieval by 10mM citrate buffer, pH6.0, 100°C for 10min.

Dilution : 1:50



GTx84860 ICC/IF Image

ICC/IF analysis of COS7 cells transiently transfected with Rho GDI alpha plasmid using GTx84860 Rho GDI alpha antibody [1A7].



GTx84860 WB Image

WB analysis of wild type and Rho GDI alpha knock out 293T cell lysate(10µg per lane) using GTx84860 Rho GDI alpha antibody [1A7].

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).