

Apolipoprotein M antibody [3H3]

Cat. No. GTX84882

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.85 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 23-188 of human APOM
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx84882 WB Image

WB analysis of various samples using GTx84882 Apolipoprotein M antibody [3H3].

Loading : 10 ug per lane


Dilution : 1:200



GTx84882 WB Image

WB analysis of HEK293T cells transfected with Apolipoprotein M plasmid (Right) or empty vector (Left) for 48 hrs using GTx84882 Apolipoprotein M antibody [3H3].

Loading : 5 ug per lane




GTx84882 WB Image

WB analysis of HepG2 cell lysate using GTx84882 Apolipoprotein M antibody [3H3].

Loading : 10 ug per lane

Dilution : 1:200



GTx84882 IHC-P Image

IHC-P analysis of liver tissue using GTx84882 Apolipoprotein M antibody [3H3].

Antigen retrieval : Heat-induced epitope retrieval by 10mM citrate buffer, pH6.0, 100°C for 10min.



For full product information, images and publications, please visit our [website](https://www.genetex.com).