

# CRISP2 antibody

**Cat. No. GTX85002**

|                     |                  |
|---------------------|------------------|
| <b>Host</b>         | Rabbit           |
| <b>Clonality</b>    | Polyclonal       |
| <b>Isotype</b>      | IgG              |
| <b>Applications</b> | WB, IHC-P, ELISA |
| <b>Reactivity</b>   | Human            |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.5 - 1 µg/mL        |
| IHC-P              | 10 µg/mL             |
| ELISA              | Assay dependent      |

Not tested in other applications.

**Calculated MW** 27 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | CRISP2 antibody was raised against a 15 amino acid synthetic peptide near the amino terminus of human CRISP2. The immunogen is located within the first 50 amino acids of CRISP2.                                          |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |


### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx85002 WB Image**

WB analysis of human testis tissue lysate using GTx85002 CRISP2 antibody.  
Working concentration : (A) 0.5 and (B) 1  $\mu$ g/ml


**GTx85002 IHC-P Image**

IHC-P analysis of human testis tissue using GTx85002 CRISP2 antibody.  
Working concentration : 20  $\mu$ g/ml


**GTx85002 IHC-P Image**

IHC-P analysis of human testis tissue using GTx85002 CRISP2 antibody.  
Working concentration : 10  $\mu$ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).