

HAAO antibody

Cat. No. GTX85015

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1µg/mL
IHC-P	2-20µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	HAAO antibody was raised against a 17 amino acid synthetic peptide near the amino terminus of human HAAO. The immunogen is located within amino acids 30 - 80 of HAAO.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx85015 IHC-P Image

IHC-P analysis of human liver tissue using GTx85015 HAAO antibody.
Working concentration : 2.5 µg/ml




GTx85015 IHC-P Image

IHC-P analysis of human liver tissue using GTx85015 HAAO antibody.
Working concentration : 20 µg/ml



GTx85015 WB Image

WB analysis of mouse liver tissue lysate using GTx85015 HAAO antibody.
Working concentration : 1 µg/ml



GTx85015 IHC-P Image

IHC-P analysis of mouse colon tissue using GTx85015 HAAO antibody.
Dilution : 2µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).